BioRem-2000

8003-005, 8003-055 BioRem-2000 Parts Cleaner Liquid™

Description

- ♦ Highly effective Microbe-Powered formula (liquid) to be used for parts' cleaning.
- Designed for use with manual, immersion, dip tank and agitation cleaning applications.
- Remediates oil and grease into carbon dioxide and water.
- Contains inhibitors for superior rust and corrosion prevention.
- Excellent rinseability, fast drying leaving no residue.
- Non-toxic, non-caustic, non-flammable and non-hazardous.
- Dilute 10:1 with water.

Technology Profile

Number of Different Microbial Strains	12
Microbial Count	50 Billion/gram
Microbial Characteristic	All GRAS Listed
Number of Enzyme Species	7
Enzyme Activity	6,000 u/mg.
pH Activity Range	5-11 pH
Appearance	Amber Liquid
Bioluminecse Test	Positive for Living Cells
Salmonella	Negative
Listeria	Negative
Phosphorous	Non-Detect

Technical Information

Usage	Dilution Ratio	RTU
	Appearance	Liquid
Physical	Color	Amber
Properties	Fragrance	None
	рН	7
	Shelf Life	Minimum 1 Year
Primary	8003-005	5 gal.
Packaging	8003-055	55 gal.



Advantages

The unique characteristic of our microbial blend is its ability to adapt to the changing distribution of hydrocarbon by-products and produce more enzymes needed to digest a particular type of chemical. The food grade botanical extracts increase the attack ratio for the microbes and enzymes.

Application

- 1. Dilute 10:1 with water. For best results heat the solution to 125 F.
- 2. Agitate for heavy, stubborn hydrocarbon build-up.
- 3. Replace when the solution becomes cloudy. (approximately 30 to 45 days).

